


X-ray absorption spectroscopy in materials science

Martin Müller

Institut für Experimentelle und Angewandte Physik
der Christian-Albrechts-Universität zu Kiel

- Photoelectric absorption
- XANES
- EXAFS
- Instrumentation
- Examples

AG M. MÜLLER
PHYSIK BIOLOGISCHER
KOMPOSITMATERIALIEN




INSTITUT FÜR EXPERIMENTELLE
UND ANGEWANDTE PHYSIK
DER CAU ZU KIEL

Uses of X-ray spectroscopy

XAFS: X-ray absorption fine structure

local technique:

- element-specific
- diluted or concentrated systems can be studied
- study of disordered systems possible
- wide temperature and pressure range



Photoelectric absorption


(a) ionisation energy of inner shell smaller than X-ray energy




photoemission
of core electron

(b,c) two possible secondary processes to fill the hole



(a) Photoelectric absorption




(b) Fluorescent X-ray emission



(c) Auger electron emission



Transmission and absorption cross section



transmitted intensity decays exponentially with thickness:

$$T = \frac{I}{I_0} = e^{-\mu z}$$

absorption
coefficient


$$\mu = \left(\frac{\rho_m N_A}{A} \right) \sigma_a$$

absorption
cross-section

mass density

atomic mass number

Avogadro's number





Absorption edges

general behaviour
of cross-section:



$$\sigma_a \propto Z^4$$

$$\sigma_a \propto E^{-3}$$

sudden increase in
absorbance at
ionisation energy
of core electron
= **absorption edge**




Absorption edges




Fine structure of absorption edges


example: krypton



2D crystalline Kr



AG M. WÖLLER
PHYSIK BIOLOGISCHER
KOMPOSITMATERIALIEN



INSTITUT FÜR EXPERIMENTELLE
UND ANGEWANDTE PHYSIK
DER CAU ZU KIE...

XAFS

X-ray **A**bsorption **N**ear-**E**dge **S**tructure

Extended **X**-ray **A**bsorption **F**ine **S**tructure

XAFS accessible elements

XANES only

EXAFS difficult

K-edge EXAFS

L3/K-edge EXAFS


L3-edge EXAFS

H																			He
Li	Be										B	C	N	O	F	Ne			
Na	Mg										Al	Si	P	S	Cl	Ar			
K	Ca	Sc	Ti	V	Cr	Mn	Fe	Co	Ni	Cu	Zn	Ga	Ge	As	Se	Br	Kr		
Rb	Sr	Y	Zr	Nb	Mo	Tc	Ru	Rh	Pd	Ag	Cd	In	Sn	Sb	Te	I	Xe		
Cs	Ba	* Lu	* Hf	* Ta	* W	* Re	* Os	* Ir	* Pt	* Au	* Hg	* Tl	* Pb	* Bi	* Po	* At	* Rn		
Fr	Ra	* Lr	* Rf	* Db	* Sg	* Bh	* Hs	* Mt	* Ds	* Rg									

* La	* Ce	* Pr	* Nd	* Pm	* Sm	* Eu	* Gd	* Tb	* Dy	* Ho	* Er	* Tm	* Yb
* Ac	* Th	* Pa	* U	* Np	* Pu	* Am	* Cm	* Bk	* Cf	* Es	* Fm	* Md	* No

X-ray absorption spectroscopy in materials science



- Photoelectric absorption
- XANES
- EXAFS
- Instrumentation
- Examples



Energy shift of edge position



energy to eject core electron depends on charge it experiences:
edge energy depends on **oxidation state**

Mn compound: shift of edge position





Pre-edge peaks

electron might also be excited into
bound states




Pre-edge peaks




E (eV)

depend on **geometry**:

- oxidation state
- site symmetry
- surrounding ligands
- nature of bonding




AG M. MÜLLER
PHYSIK BIOLOGISCHER
KOMPOSITMATERIALIEN




INSTITUT FÜR EXPERIMENTELLE
UND ANGEWANDTE PHYSIK
DER SAU ZU KIE...

XANES white lines



in transition metals:
area of white line
indicates number of
empty d-states




XANES fingerprinting

XANES characteristic of *chemical environment* and *valence state*:


- fingerprinting
- phase analysis by **linear combination** of known species

unknown
S spectrum

known references
(oxidation numbers
0, 2, 4, 6)




high sensitivity
to electronic states!




X-ray absorption spectroscopy in materials science

- Photoelectric absorption
- XANES
- EXAFS
- Instrumentation
- Examples




Principle of EXAFS oscillations


Interference of wavefunctions of photoelectron and of **backscattered** electrons (from neighbouring atoms)



AG M. WOLLER
PHYSIK BIOLOGISCHER
KOMPOSITMATERIALIEN





INSTITUT FÜR EXPERIMENTELLE
UND ANGEWANDTE PHYSIK
DER SAU ZU KIEL



Animation

Parameters accessible with EXAFS

- **type** of atoms surrounding central absorber ($Z \pm 3$)
- **number** of atoms surrounding absorber ($\pm 20\%$)
- **distances** absorber – scatterer (accuracy 0.1 Å)



EXAFS data analysis

EXAFS signal:

$$\chi(q(E)) = \frac{\mu_{\chi}(E) - \mu_0(E)}{\mu_0(E)}$$


isolated atom

in material of interest



electron wave vector

$$\frac{\hbar^2 q^2}{2m} = E - \hbar\omega_K$$

Fourier transform:
radial distribution function



CdTe nanocrystals





The EXAFS formula ...

$$\chi(q) \propto \sum_j N_j \frac{t_j(q) \sin(2qR_j + \delta_j(q))}{R_j^2} e^{-2(q\sigma_j)^2} e^{-2R_j/\Lambda}$$


- sum over j **neighbouring** shells
 - **goal**: extract radii R_j and occupation numbers N_j
 - **damping** due to loss (mean free path Λ) and disorder (Debye-Waller factor)
 - phase shifts
 - backscattering amplitude
- } difficult...

CdTe bulk





X-ray absorption spectroscopy in materials science

- Photoelectric absorption
- XANES
- EXAFS
- Instrumentation
- Examples



Absorption spectrometer



(fast scanning with piezos: QEXAFS)



IC: ionisation chamber


Fluorescence

scanning




EXAFS / XANES in one shot: DEXAFS

„white“
synchrotron
radiation




polychromator:
bent Si(111) crystal

energy dispersion
translated into position
dependence on detector




energy-dispersive

AG M. WOLLER
PHYSIK BIOLOGISCHER
KOMPOSITMATERIALIEN




INSTITUT FÜR EXPERIMENTELLE
UND ANGEWANDTE PHYSIK
DER CAU ZU KASSEL

Bent crystal polychromator @ ID24 (ESRF)



AG M. MÖLLER
PHYSIK BIOLOGISCHER
KOMPOSITMATERIALIEN




INSTITUT FÜR EXPERIMENTELLE
UND ANGEWANDTE PHYSIK
DER CAU ZU KIE...


Alternatives to transmission measurements

fluorescence
(high sensitivity)


Auger yield
(if fluorescence yield low)




(a) Photoelectric absorption



(b) Fluorescent X-ray emission




(c) Auger electron emission






X-ray absorption spectroscopy in materials science

- Photoelectric absorption
- XANES
- EXAFS
- Instrumentation
- Examples




Homogeneous catalytic reaction mechanisms

Cu(II) catalyst (arylation)



- **XANES** time resolution: 10-200 ms
- combination with UV-vis spectroscopy

AG M. WÖLLER
PHYSIK BIOLOGISCHER
KOMPOSITMATERIALIEN




INSTITUT FÜR EXPERIMENTELLE
UND ANGEWANDTE PHYSIK
DER CAU ZU KIE...

Adsorbate-induced phase change in Rh catalysts

car exhaust **catalytic converters**:
NO converted to N₂


time-resolved **EXAFS** study



Rh **metal**


Rh nanoparticles **rapidly**
change upon exposure to
NO

Rh₂O₃ **oxide**




Synchronizing IR spectroscopy and XAFS

again de-NOx Rh catalyst: many different components



sub-second **EXAFS** and IR time resolution

AG M. WÖLLER
PHYSIK BIOLOGISCHER
KOMPOSITMATERIALIEN



INSTITUT FÜR EXPERIMENTELLE
UND ANGEWANDTE PHYSIK
DER CAU ZU KIEL